

## Product datasheet for **MP213966**

### Rptor Mouse qPCR Primer Pair (NM\_028898)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	74370
Forward Sequence:	CTTCCTATCCGTCTTGCCAGAC
Reverse Sequence:	CTCCAGACAGATGGCAATCAGG
ACCN:	<a href="#">NM_028898</a> , <a href="#">NM_028898.1</a> , <a href="#">NM_028898.2</a> , <a href="#">NM_028898.3</a> , <a href="#">BC023733</a> , <a href="#">BC052202</a> , <a href="#">BC059877</a> , <a href="#">BC060140</a> , <a href="#">BC118971</a> , <a href="#">BC148309</a> , <a href="#">BC167200</a>
Synonyms:	4932417H02Rik; mKIAA1303; r; Rap; Raptor
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)