

Product datasheet for **MP213977**

TxIng Mouse qPCR Primer Pair (NM_178935)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	353170
Forward Sequence:	CCACAATGCCAAACTGCGACAG
Reverse Sequence:	GCTGAAGTTTGGCATCCACAAGC
ACCN:	<u>BC117898</u> , <u>BC117899</u> , <u>NM_178935</u> , <u>NM_178935.1</u> , <u>NM_178935.2</u> , <u>NM_178935.3</u> , <u>NM_178935.4</u> , <u>NM_178935.5</u> , <u>BC059099</u> , <u>BC080778</u> , <u>BQ128960</u> , <u>BY602007</u>
UniProt ID:	<u>Q8BHN1</u>
Synonyms:	4932441K18Rik; Fiat; Lrp
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)