

Product datasheet for **MP214084**

Arhgap27 Mouse qPCR Primer Pair (NM_183288)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	544817
Forward Sequence:	CGCCTAGCGACTCGGAGAATG
Reverse Sequence:	ACTCTGCTGGACTGTGACCTCC
ACCN:	NM_183288 , NM_183288.3 , BC138292 , BC005635 , BC096671 , BC138293 , BY552449
UniProt ID:	A2AB59
Synonyms:	2310069I04Rik; 5730442P18Rik; Gm21552; Sh3d20
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)