

## Product datasheet for **MP214711**

### Cgas Mouse qPCR Primer Pair (NM\_173386)

#### Product data:

|                   |   |
|-------------------|---|
| Product Type:     | qPCR Primer Pairs   |
| Gene ID:          | 214763  |
| Forward Sequence: | TTCCACGAGGAAATCCGCTGAG  |
| Reverse Sequence: | CAGCAGGGCTTCCTGGTTTTTC  |
| ACCN:             | <u><a href="#">BC052196</a></u> , <u><a href="#">NM_173386</a></u> , <u><a href="#">NM_173386.1</a></u> , <u><a href="#">NM_173386.3</a></u> , <u><a href="#">NM_173386.4</a></u> , <u><a href="#">NM_173386.5</a></u> , <u><a href="#">BC005787</a></u> ,<br><u><a href="#">BC044854</a></u> , <u><a href="#">BC145651</a></u> , <u><a href="#">BC145653</a></u> |
| UniProt ID:       | <u><a href="#">Q8C6L5</a></u>   |
| Synonyms:         | E330016A19Rik; Mb21d1   |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)  |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.                      |
| Storage:          | Store at -20°C.   |
| Stability:        | The primer mix is stable for one year from date of shipping.  |



[View online »](#)