

Product datasheet for **MP214777**

Rnf14 Mouse qPCR Primer Pair (NM_020012)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	56736
Forward Sequence:	TGCTTTCTGCCACCACTTGTGC
Reverse Sequence:	GGTTGTCTAGGTGCTTGCAGAG
ACCN:	<u>BC054841</u> , <u>BC094250</u> , <u>NM_020012</u> , <u>NM_020012.1</u> , <u>NM_020012.2</u> , <u>BC049079</u> , <u>BC052205</u>
UniProt ID:	<u>Q9JI90</u>
Synonyms:	2310075C09Rik; 2610005D23Rik; AA986456; AU041447; D7Bwg0165e; D18Ert188e; Triad2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)