

Product datasheet for **MP214778**

Rnf141 Mouse qPCR Primer Pair (NM_025999)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67150
Forward Sequence:	GCAAAGCACGTCACCTTGGTTC
Reverse Sequence:	TTCCAGAAGGCAGAGGAGTCAG
ACCN:	<u>BC018553</u> , <u>NM_025999</u> , <u>NM_025999.1</u> , <u>NM_025999.2</u> , <u>NM_025999.3</u> , <u>BC051516</u> , <u>BC096563</u> , <u>BY265291</u>
UniProt ID:	<u>Q99MB7</u>
Synonyms:	2610110L04Rik; AA792898; AU022812; ZFP36; ZNF230
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)