

Product datasheet for **MP214783**

Rnfl49 Mouse qPCR Primer Pair (NM_001033135)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67702
Forward Sequence:	GAGAGGAGTGAAGCAACCTTCC
Reverse Sequence:	TCCTCCATGCTGAGGCTCACTC
ACCN:	<u>BC115968</u> , <u>NM_001033135</u> , <u>NM_001033135.1</u> , <u>NM_001033135.2</u> , <u>NM_001033135.3</u> , <u>BE627920</u> , <u>BF018252</u> , <u>BY743657</u>
UniProt ID:	<u>Q3U2C5</u>
Synonyms:	1600023E10Rik; AV037261; C76763; Gm15832; Greul4; OTTMUSG00000026591
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

