

## Product datasheet for **MP214798**

### Rnf2 Mouse qPCR Primer Pair (NM\_011277)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	19821
Forward Sequence:	CGCATCAGGAAAGGGTCTTAGC
Reverse Sequence:	CTATCTGCTGCTTTTTGCCTCGC
ACCN:	<u><a href="#">BC020122</a></u> , <u><a href="#">NM_011277</a></u> , <u><a href="#">NM_011277.1</a></u> , <u><a href="#">NM_011277.2</a></u> , <u><a href="#">NM_011277.3</a></u>
UniProt ID:	<u><a href="#">Q9CQJ4</a></u>
Synonyms:	AI326319; AI450156; AU019207; dinG; Ring1B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)