

## Product datasheet for **MP214945**

### S100a10 Mouse qPCR Primer Pair (NM\_009112)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20194
Forward Sequence:	GACAAAGGAGGACCTGAGAGTG
Reverse Sequence:	CTCTGGAAGCCCACTTTGCCAT
ACCN:	<a href="#">BC025044</a> , <a href="#">NM_009112</a> , <a href="#">NM_009112.1</a> , <a href="#">NM_009112.2</a> , <a href="#">BU559648</a>
UniProt ID:	<a href="#">P08207</a>
Synonyms:	42C; AA409961; AL024248; Cal1l; CAL12; CLP11; p10; p11
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)