

Product datasheet for **MP214956**

Calprotectin (S100a9) Mouse qPCR Primer Pair (NM_009114)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20202
Forward Sequence:	TGGTGAAGCACAGTTGGCAAC
Reverse Sequence:	CAGCATCATACACTCCTCAAAGC
ACCN:	<u>BC027635, NM_009114, NM_009114.1, NM_009114.2, NM_009114.3, BF732133, BY112406</u>
UniProt ID:	<u>P31725</u>
Synonyms:	60B8Ag; AW546964; BEE22; Cagb; GAGB; L1Ag; MRP14; p14
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)