

Product datasheet for **MP214968**

Amd1 Mouse qPCR Primer Pair (NM_009665)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	11702
Forward Sequence:	GCAAAGGATGTCACTCGTGAGAG
Reverse Sequence:	GTTCTGGAGTGATGTGAATAGTCC
ACCN:	BC011110 , BC071220 , BC080791 , BC092072 , NM_009665 , NM_009665.1 , NM_009665.2 , NM_009665.3 , NM_009665.4 , NM_009665.5 , BC086902 , BC138696 , BC138697
UniProt ID:	P0DMN7
Synonyms:	AdoMetDC; adoMetDC1; Amd-1; SAMDC; SAMDC 1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)