

## Product datasheet for **MP214989**

### Sar1b Mouse qPCR Primer Pair (NM\_025535)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66397
Forward Sequence:	GGATTGTGCAGACCATGAAAGGC
Reverse Sequence:	CACTGATGGCTTCAGGTCTGTC
ACCN:	<u><a href="#">BC082550</a></u> , <u><a href="#">BC085178</a></u> , <u><a href="#">NM_025535</a></u> , <u><a href="#">NM_025535.1</a></u> , <u><a href="#">NM_025535.2</a></u> , <u><a href="#">BB749343</a></u> , <u><a href="#">BY068863</a></u>
UniProt ID:	<u><a href="#">Q9CQC9</a></u>
Synonyms:	2310075M17Rik; 2900019I22Rik; CMRD; Sara1b; Sara2; Sarb
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)