

Product datasheet for **MP215005**

Scarb1 Mouse qPCR Primer Pair (NM_016741)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20778
Forward Sequence:	ACACCCGAATCCTCGCTGGAAT
Reverse Sequence:	CCGTTGGCAAACAGAGTATCGG
ACCN:	<u>BC004656</u> , <u>NM_016741</u> , <u>NM_016741.1</u> , <u>NM_016741.2</u> , <u>BQ176029</u> , <u>BY018765</u>
UniProt ID:	<u>Q61009</u>
Synonyms:	A1120173; CD36; Cd36l1; Chohd1; Cla-1; Cla1; D5Ertd460e; Hdlq1; Hlb398; mSR-BI; SR-B1; SR-BI
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)