

Product datasheet for **MP215070**

Secisbp2 Mouse qPCR Primer Pair (NM_029279)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	75420
Forward Sequence:	CCGAGTTATGAAGTCCTGACGG
Reverse Sequence:	CCTCTGTGTCTTCTTCCGACG
ACCN:	<u>BC125391</u> , <u>BC125393</u> , <u>NM_029279</u> , <u>NM_029279.1</u> , <u>NM_029279.2</u> , <u>BC038877</u> , <u>BC144795</u> , <u>BP760396</u>
Synonyms:	2210413N07Rik; 2810012K13Rik; SB; SBP2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)