

## Product datasheet for **MP215080**

### Spp1 Mouse qPCR Primer Pair (NM\_009263)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20750
Forward Sequence:	GCTTGGCTTATGGACTGAGGTC
Reverse Sequence:	CCTTAGACTCACCGCTCTTCATG
ACCN:	<a href="#">BC002113</a> , <a href="#">BC057858</a> , <a href="#">BC080720</a> , <a href="#">NM_009263</a> , <a href="#">NM_009263.1</a> , <a href="#">NM_009263.2</a> , <a href="#">NM_009263.3</a> , <a href="#">BC014284</a> , <a href="#">BC020355</a> , <a href="#">BC029991</a>
UniProt ID:	<a href="#">P10923</a>
Synonyms:	2AR; Apl-1; BNSP; Bsp; BSPI; Eta; ETA-1; OP; Opn; Opnl; Ric; Spp-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)