

Product datasheet for **MP215184**

Serpina3g Mouse qPCR Primer Pair (NM_009251)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20715
Forward Sequence:	GCATCAGGGAAGTCTTCTCCAC
Reverse Sequence:	CACAACATCCGACACCTGCCAT
ACCN:	BC119537 , BC119538 , NM_009251 , NM_009251.1 , BC002065 , BC057144 , BC096664 , BE850123 , BE851987 , NM_009251.2
UniProt ID:	Q5I2A0
Synonyms:	2A2; A1119734; spi2; Spi2-1; Spi2/eb.1; Spi2A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)