

Product datasheet for **MP215262**

Stk11 Mouse qPCR Primer Pair (NM_011492)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20869
Forward Sequence:	GCCTGGAATACCTACACAGCCA
Reverse Sequence:	GCAGGTGTCATCCACAGCGAAA
ACCN:	BC052379 , NM_011492 , NM_011492.1 , NM_011492.2 , NM_011492.3 , NM_011492.4 , BE951985 , BF468348 , NM_011492.5
UniProt ID:	Q9WTK7
Synonyms:	AA408040; Lkb1; Par-4; R75140
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)