

Product datasheet for **MP215285**

Nemf Mouse qPCR Primer Pair (NM_025441)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66244
Forward Sequence:	GCTGCACAAGTGAAGTCATGGC
Reverse Sequence:	GAGTTCTACGTCTCCAGGCATTC
ACCN:	<u>NM_025441</u> , <u>NM_025441.3</u> , <u>BB337197</u> , <u>BC027272</u> , <u>BC037106</u> , <u>BC039647</u> , <u>BC053488</u> , <u>BC080716</u> , <u>BF452703</u> , <u>BG082920</u> , <u>BY241581</u> , <u>BY741921</u>
UniProt ID:	<u>Q8CCP0</u>
Synonyms:	1500011112Rik; 4933405E14Rik; Sdccag1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)