

## Product datasheet for **MP215317**

### Setd7 Mouse qPCR Primer Pair (NM\_080793)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	73251
Forward Sequence:	TTGACGGAGAGATGCTCGAAGG
Reverse Sequence:	GAAGGAGAGCATCGCTGGAGAT
ACCN:	<a href="#">BC050190</a> , <a href="#">NM_080793</a> , <a href="#">NM_080793.1</a> , <a href="#">NM_080793.2</a> , <a href="#">NM_080793.3</a> , <a href="#">NM_080793.4</a> , <a href="#">NM_080793.5</a> , <a href="#">BC008522</a>
UniProt ID:	<a href="#">Q8VHL1</a>
Synonyms:	1600028F23Rik; H3K4MT; KMT7; mKIAA1717; Set7; Set7/9
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)