

Product datasheet for **MP215438**

Stat5a Mouse qPCR Primer Pair (NM_011488)

Product data:

| | |
|-------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 20850 |
| Forward Sequence: | CCTGTTTGAGTCTCAGTTCAGCG |
| Reverse Sequence: | TGGCAGTAGCATTGTGGTCCTG |
| ACCN: | BC008998 , BC013274 , BC056647 , NM_011488 , NM_011488.1 , NM_011488.2 , NM_011488.3 , BY277646 |
| UniProt ID: | P42230 |
| Synonyms: | AA959963; STAT5 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)