

Product datasheet for **MP215512**

Steap2 Mouse qPCR Primer Pair (NM_028734)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	74051
Forward Sequence:	GAGCAACGCTTTGAACTGGAGAG
Reverse Sequence:	GGCAAGAACGAAGTTTGGTGGTG
ACCN:	<u>NM_028734</u> , <u>NM_028734.2</u> , <u>NM_028734.3</u> , <u>NM_028734.4</u> , <u>NM_028734.5</u> , <u>BC150997</u> , <u>BB751127</u> , <u>BC043459</u> , <u>BC150881</u> , <u>BC150889</u> , <u>BC151008</u>
UniProt ID:	<u>Q8BWB6</u>
Synonyms:	4921538B17Rik; AI930049; AW045895; IPCA-1; IPCA1; PCANAP1; STAMP1; STMP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)