

## Product datasheet for **MP215553**

### Sgta Mouse qPCR Primer Pair (NM\_024499)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	52551
Forward Sequence:	TGGAGGACAGTGACCTTGCTCT
Reverse Sequence:	TCTTCAGAGGGTGGTGCCTGT
ACCN:	<a href="#">BC003836</a> , <a href="#">NM_024499</a> , <a href="#">NM_024499.1</a> , <a href="#">NM_024499.2</a>
UniProt ID:	<a href="#">Q8BJU0</a>
Synonyms:	5330427H01Rik; AI194281; D10ErtD190e; Sgt; Stg
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)