

Product datasheet for **MP215598**

Smo Mouse qPCR Primer Pair (NM_176996)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	319757
Forward Sequence:	GAGGCTACTTCCTCATCAGAGG
Reverse Sequence:	GCTGAAGGTGATGAGCACAAAGC
ACCN:	BC048091 , BC096028 , NM_176996 , NM_176996.1 , NM_176996.2 , NM_176996.3 , NM_176996.4 , BB659739 , BB665395 , BC030377 , BQ947565
UniProt ID:	P56726
Synonyms:	bnb; E130215L21Rik; Smoh; smoothed
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)