

Product datasheet for **MP215622**

Scn4b Mouse qPCR Primer Pair (NM_001013390)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	399548
Forward Sequence:	GGCAGATACACCTGCTTCGTGA
Reverse Sequence:	TGAGAGTCACCGTGTGTCCAC
ACCN:	<u>NM_001013390</u> , <u>NM_001013390.1</u> , <u>NM_001013390.2</u> , <u>BC028468</u> , <u>BC050045</u> , <u>BC139399</u> , <u>BC139400</u> , <u>BK001031</u> , <u>NM_001013390.3</u>
UniProt ID:	<u>Q7M729</u>
Synonyms:	Gm1471
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)