

Product datasheet for **MP215637**

Slc1a3 Mouse qPCR Primer Pair (NM_148938)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20512
Forward Sequence:	GCGATTGGTCGCGGTGATAATG
Reverse Sequence:	CGACAATGACTGTCACGGGTAC
ACCN:	<u>BC058711</u> , <u>BC066154</u> , <u>NM_148938</u> , <u>NM_148938.1</u> , <u>NM_148938.2</u> , <u>NM_148938.3</u> , <u>BC028638</u> , <u>BC055300</u> , <u>BY281420</u>
UniProt ID:	<u>P56564</u>
Synonyms:	AI504299; B430115D02Rik; Eaat1; GLAST; GLAST-1; GLU-T; GluT-1; Gmt1; MGluT1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)