

Product datasheet for **MP215659**

Slc12a7 Mouse qPCR Primer Pair (NM_011390)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20499
Forward Sequence:	AACTGGCGTCCACGCTTCAAGT
Reverse Sequence:	GGCAATGAGCATGGCGAAAAGG
ACCN:	<u>NM_011390</u> , <u>NM_011390.1</u> , <u>NM_011390.2</u> , <u>BC141107</u> , <u>BC051558</u> , <u>BC053470</u> , <u>BC059242</u> , <u>BC060133</u> , <u>BC079677</u>
UniProt ID:	<u>Q9WVL3</u>
Synonyms:	AA408796; D13Ertd261e; Kcc4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)