

Product datasheet for **MP215694**

Slc18a2 Mouse qPCR Primer Pair (NM_172523)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	214084
Forward Sequence:	CCTCTTACGACCTTGCTGAAGG
Reverse Sequence:	GCTGCCACTTTCGGGAACACAT
ACCN:	<u>BC078449</u> , <u>NM_172523</u> , <u>NM_172523.1</u> , <u>NM_172523.2</u> , <u>NM_172523.3</u>
UniProt ID:	<u>Q8BRU6</u>
Synonyms:	1110037L13Rik; 9330105E13; Vmat2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)