

Product datasheet for **MP215723**

Slc22a18 Mouse qPCR Primer Pair (NM_008767)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	18400
Forward Sequence:	GTCCTCAGCTTTACCTGTGTCC
Reverse Sequence:	CTTGACCAGGAATACAGGCAGC
ACCN:	BC003734 , BC015078 , NM_008767 , NM_008767.1 , NM_008767.2 , BB869719 , BE865039 , BU758157 , NM_008767.3
UniProt ID:	Q78KK3
Synonyms:	AW260131; BWR1A; BWSCR1A; HET; Impt1; ITM; Orctl2; p45-BWR1A; Slc22a1l; TSSC5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)