

Product datasheet for **MP215755**

Slc25a25 Mouse qPCR Primer Pair (NM_146118)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	227731
Forward Sequence:	GGAACAGCAGGCGGAGAAGATT
Reverse Sequence:	CACTCGTTCAGTCGATGGTCA
ACCN:	<u>BC019978</u> , <u>BC022114</u> , <u>BC037109</u> , <u>BC066998</u> , <u>NM_146118</u> , <u>NM_146118.1</u> , <u>NM_146118.2</u> , <u>NM_146118.3</u> , <u>BC019978.1</u> , <u>BB653459</u> , <u>BC004720</u> , <u>BY772900</u>
UniProt ID:	<u>A2ASZ8</u>
Synonyms:	1110030N17Rik; MCSC; mKIAA1896
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)