

Product datasheet for **MP215822**

Slc35a2 Mouse qPCR Primer Pair (NM_078484)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	22232
Forward Sequence:	CCTCCAGTATGTTGCCATCAGC
Reverse Sequence:	GCTGCGATTCAACATGAGCACAG
ACCN:	<u>BC037701</u> , <u>NM_078484</u> , <u>NM_078484.1</u> , <u>NM_078484.2</u> , <u>BX512781</u>
UniProt ID:	<u>Q9R0M8</u>
Synonyms:	AI327289; Had-1; Had1; Sfc8; Ugalt; UGT
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)