

## Product datasheet for **MP215832**

### Slc35c2 Mouse qPCR Primer Pair (NM\_144893)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	228875
Forward Sequence:	CTGCTTCTCCATAGGCATCACG
Reverse Sequence:	GCCTTGTGGCTGGAGCATTGAA
ACCN:	<a href="#">BC018327</a> , <a href="#">BC094025</a> , <a href="#">NM_144893</a> , <a href="#">NM_144893.1</a> , <a href="#">NM_144893.2</a> , <a href="#">BB772752</a> , <a href="#">BB849961</a> , <a href="#">BC064092</a> , <a href="#">BE952356</a>
UniProt ID:	<a href="#">Q8VCX2</a>
Synonyms:	C85957; CGI-15; D2Wsu58e; Ovcov1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)