

Product datasheet for **MP215892**

Slc44a2 Mouse qPCR Primer Pair (NM_152808)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	68682
Forward Sequence:	TGCTCTACCCACTGGTGACCTT
Reverse Sequence:	GCAGGTTTTCTCAGAAGTGGG
ACCN:	BC031535 , NM_152808 , NM_152808.1 , NM_152808.2 , NM_152808.3 , BC013213 , BC023763 , BU756489 , BY339157
UniProt ID:	Q8BY89
Synonyms:	1110028E10Rik; CTL2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)