

Product datasheet for **MP215908**

Slc5a4b Mouse qPCR Primer Pair (NM_023219)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64454
Forward Sequence:	CCTCTTCACCATCGACCTGTAC
Reverse Sequence:	CACAGGATGCTGACGGCGATTA
ACCN:	<u>BC127615</u> , <u>NM_023219</u> , <u>NM_023219.1</u> , <u>NM_023219.2</u> , <u>BC109133</u> , <u>BC109134</u> , <u>BC127612</u> , <u>BC128465</u>
Synonyms:	2010104G07Rik; pSGLT2; SAAT1; SGLT3-b; SGLT3b
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)