

Product datasheet for **MP215957**

Slc9a2 Mouse qPCR Primer Pair (NM_001033289)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	226999
Forward Sequence:	GTGTCCACTGTTGGGAAGAACC
Reverse Sequence:	ACTGGTCCTGAAGGTTAGCGG
ACCN:	<u>BC103685</u> , <u>BC104737</u> , <u>NM_001033289</u> , <u>NM_001033289.1</u> , <u>NM_001033289.2</u> , <u>BC105676</u>
Synonyms:	2210416H12Rik; 4932415O19Rik; AV333884; NHE2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)