

## Product datasheet for **MP216041**

### Spg7 Mouse qPCR Primer Pair (NM\_153176)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	234847
Forward Sequence:	GGTTCTCCAACACGGAAGAGGA
Reverse Sequence:	GCCTCATCAGAGCATTGTCCAG
ACCN:	<u><a href="#">NM_153176</a></u> , <u><a href="#">NM_153176.1</a></u> , <u><a href="#">NM_153176.2</a></u> , <u><a href="#">NM_153176.3</a></u> , <u><a href="#">NM_153176.4</a></u> , <u><a href="#">BC138141</a></u> , <u><a href="#">BC024466</a></u> , <u><a href="#">BC024986</a></u> , <u><a href="#">BC051051</a></u> , <u><a href="#">BC055488</a></u> , <u><a href="#">BC096690</a></u>
UniProt ID:	<u><a href="#">Q3ULF4</a></u>
Synonyms:	AI452278; AU015315; Cmar; PGN
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)