

Product datasheet for **MP216053**

Sptbn1 Mouse qPCR Primer Pair (NM_175836)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20742
Forward Sequence:	GTGGGAGAACAGGCAAAACCTC
Reverse Sequence:	GGTGGCATTTCAGTATGAGCC
ACCN:	NM_175836 , NM_175836.1 , NM_175836.2 , BC023161 , BC024833 , BC051653 , BC059940 , BC092544 , BC156448 , BC172699
UniProt ID:	Q62261
Synonyms:	9930031C03Rik; AL033301; elf1; elf3; mKIAA4049; Spnb-2; Spnb2; SPTB2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)