

Product datasheet for **MP216074**

Spag9 Mouse qPCR Primer Pair (NM_027569)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	70834
Forward Sequence:	TGCTGGTTTGGAACTGAAGGC
Reverse Sequence:	CCAGACCTGACTGGACAACCTCT
ACCN:	BC060100 , BC094670 , NM_027569 , NM_027569.1 , NM_027569.2 , BC094670.1 , BB618351 , BC019555 , BC052058 , BC060506 , BC096410 , NM_027569.3
UniProt ID:	Q58A65
Synonyms:	3110018C07Rik; 4733401I23Rik; 4831406C20Rik; AW552012; Jip4; JLP; JSAP2; JSAP2a
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)