

Product datasheet for **MP216099**

Spata6 Mouse qPCR Primer Pair (NM_026470)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 67946 |
| Forward Sequence: | GTGGTGCTGAAAGACAAGGAGG |
| Reverse Sequence: | CACATCTCCAGGGTCTACTGCT |
| ACCN: | <u>NM_026470</u> , <u>NM_026470.1</u> , <u>NM_026470.2</u> , <u>NM_026470.3</u> , <u>BC028429</u> , <u>BC040024</u> , <u>BC141421</u> , <u>BC146567</u> |
| UniProt ID: | <u>Q3U6K5</u> |
| Synonyms: | 1700062C23Rik; AI790763; Hash; KRP; Mash |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)