

Product datasheet for **MP216156**

Sf3b2 Mouse qPCR Primer Pair (NM_030109)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 319322 |
| Forward Sequence: | GAGCACAAGGACAGTGACGATG |
| Reverse Sequence: | GCTGCTTGAGTTCTGCCACAGT |
| ACCN: | <u>BC049118</u> , <u>NM_030109</u> , <u>NM_030109.1</u> , <u>NM_030109.2</u> , <u>BC026365</u> , <u>BC029745</u> , <u>BC045207</u> , <u>BI134679</u> , <u>BY785284</u> |
| Synonyms: | 145kDa; 2610311M13Rik; 2810441F20Rik; B230398H18Rik; SAP145; SF3b1; SF3b145; SF3b150 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)