

## Product datasheet for **MP216163**

### Srsf1 Mouse qPCR Primer Pair (NM\_173374)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	110809
Forward Sequence:	GTTTACCGAGATGGCACTGGTG
Reverse Sequence:	ACCCGGATGTAGGCAGTTTCTC
ACCN:	<u><a href="#">BC046773</a></u> , <u><a href="#">BC058627</a></u> , <u><a href="#">NM_173374</a></u> , <u><a href="#">NM_173374.1</a></u> , <u><a href="#">NM_173374.2</a></u> , <u><a href="#">NM_173374.3</a></u> , <u><a href="#">NM_173374.4</a></u> , <u><a href="#">BC051213</a></u> , <u><a href="#">BY743476</a></u>
UniProt ID:	<u><a href="#">Q6PDM2</a></u>
Synonyms:	1110054N12Rik; 5730507C05Rik; 6330415C05Rik; AI482334; Asf; AW491331; Sf; Sf2; Sfrs1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)