

Product datasheet for **MP216176**

Srsf6 Mouse qPCR Primer Pair (NM_026499)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	67996
Forward Sequence:	AGAAGGTCTCGGAGTAGGAGTC
Reverse Sequence:	GGATTTCTGCCTTTTGAGCGG
ACCN:	BC012039 , NM_026499 , NM_026499.1 , NM_026499.2 , NM_026499.3 , NM_026499.4 , BC029312 , BY783855
UniProt ID:	Q3TWW8
Synonyms:	1210001E11Rik; AI314910; AW146126; Sf; Sfrs6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)