

Product datasheet for **MP216211**

Sox11 Mouse qPCR Primer Pair (NM_009234)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20666
Forward Sequence:	GGACCTGGATTCCTTCAGTGAG
Reverse Sequence:	GTGAACACCAGGTCGGAGAAGT
ACCN:	<u>BC062095</u> , <u>BC078643</u> , <u>NM_009234</u> , <u>NM_009234.1</u> , <u>NM_009234.2</u> , <u>NM_009234.3</u> , <u>NM_009234.4</u> , <u>NM_009234.5</u> , <u>NM_009234.6</u> , <u>BC053503</u> , <u>BK001484</u>
UniProt ID:	<u>Q7M6Y2</u>
Synonyms:	1110038H03Rik; 6230403H02Rik; AI836553; end1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)