

Product datasheet for **MP216213**

Sox13 Mouse qPCR Primer Pair (NM_011439)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20668
Forward Sequence:	CCTATTCAGCCCATTCCCTGCA
Reverse Sequence:	CTTGGCTGTGAGGTTCA GTGGT
ACCN:	<u>BC007130</u> , <u>NM_011439</u> , <u>NM_011439.1</u> , <u>NM_011439.2</u> , <u>BE949927</u> , <u>BE984765</u>
UniProt ID:	<u>Q04891</u>
Synonyms:	AA407916; AW259412; AW540933; Sox-13
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)