

Product datasheet for **MP216231**

St3gal4 Mouse qPCR Primer Pair (NM_009178)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20443
Forward Sequence:	TGCCTCCAACAAGAAGCAGACC
Reverse Sequence:	AGCATCCGCTTGATGGCGATAG
ACCN:	<u>BC011121</u> , <u>BC050773</u> , <u>NM_009178</u> , <u>NM_009178.1</u> , <u>NM_009178.2</u> , <u>NM_009178.3</u> , <u>NM_009178.4</u>
UniProt ID:	<u>Q91Y74</u>
Synonyms:	SIAT4-C; Siat4c; ST3Gal IV; ST3GalIV
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)