

Product datasheet for **MP216270**

Slk Mouse qPCR Primer Pair (NM_009289)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20874
Forward Sequence:	ACGGAGGAAGTGCCTTTCTCA
Reverse Sequence:	TCAGCATCTCTGCTCAAAGCG
ACCN:	<u>BC131674</u> , <u>BC131675</u> , <u>NM_009289</u> , <u>NM_009289.1</u> , <u>NM_009289.2</u> , <u>NM_009289.3</u> , <u>BC025855</u> , <u>BC048755</u>
UniProt ID:	<u>O54988</u>
Synonyms:	9A2; AV021402; AW411554; Etk4; mKIAA0204; mSLK; SMAK; Stk2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)