

Product datasheet for **MP216303**

Srebf2 Mouse qPCR Primer Pair (NM_033218)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20788
Forward Sequence:	AGAAAGAGCGGTGGAGTCCTTG
Reverse Sequence:	GAAGTCTGGAGAATGGTGAGG
ACCN:	<u>NM_033218</u> , <u>NM_033218.1</u> , <u>BC005777</u> , <u>BC014730</u> , <u>BC023333</u> , <u>BC026556</u> , <u>BC038284</u> , <u>BC048071</u> , <u>BC069989</u> , <u>BC167259</u> , <u>BG800576</u> , <u>BQ963055</u> , <u>BQ965109</u>
UniProt ID:	<u>Q3U1N2</u>
Synonyms:	AI608257; bHLHd2; lop13; nuc; SREBP-2; SREBP2; SREBP2gc
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)