

## Product datasheet for **MP216309**

### Stra8 Mouse qPCR Primer Pair (NM\_009292)

#### Product data:

|                   |   |
|-------------------|---|
| Product Type:     | qPCR Primer Pairs   |
| Gene ID:          | 20899   |
| Forward Sequence: | CAGCTTAGAGGAGGTCAAGGAAG   |
| Reverse Sequence: | AGCATCTGGTCCAACAGCCTCA  |
| ACCN:             | <u><a href="#">BC103590</a></u> , <u><a href="#">BC103591</a></u> , <u><a href="#">NM_009292</a></u> , <u><a href="#">NM_009292.1</a></u> , <u><a href="#">BC100746</a></u> , <u><a href="#">BC103603</a></u> , <u><a href="#">NM_009292.2</a></u>  |
| UniProt ID:       | <u><a href="#">P70278</a></u>   |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec. |
| Storage:          | Store at -20 °C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.  |



[View online »](#)