

Product datasheet for **MP216350**

Sdha Mouse qPCR Primer Pair (NM_023281)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66945
Forward Sequence:	GAGATACGCACCTGTTGCCAAG
Reverse Sequence:	GGTAGACGTGATCTTTCTCAGGG
ACCN:	BC031849 , NM_023281 , BC011301
UniProt ID:	Q8K2B3
Synonyms:	1500032O14Rik; 2310034D06Rik; 4921513A11; C81073; FP; SDH2; SDHF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)