

Product datasheet for **MP216356**

Sucla2 Mouse qPCR Primer Pair (NM_011506)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 20916 |
| Forward Sequence: | GGTGTCTCTGTTCCCAAAGGCT |
| Reverse Sequence: | TTTCCTCTGCCGCGCCAAAA |
| ACCN: | BC052774 , BC056353 , BC057605 , NM_011506 , NM_011506.1 , NM_011506.2 , NM_011506.3 , BC052774.1 , BC038497 , NM_011506.4 |
| UniProt ID: | Q9Z2I9 |
| Synonyms: | 4930547K18Rik; A-SCS |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)