

## Product datasheet for **MP216357**

### Suclg2 Mouse qPCR Primer Pair (NM\_011507)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	20917
Forward Sequence:	AGCTCAAGGTGCCACTGGTAGT
Reverse Sequence:	GCTTTCTTGGCTGCATCCTCCA
ACCN:	<a href="#">BC080781</a> , <a href="#">NM_011507</a> , <a href="#">NM_011507.1</a> , <a href="#">NM_011507.2</a> , <a href="#">NM_011507.3</a> , <a href="#">BC030947</a> , <a href="#">BC043312</a> , <a href="#">BC054425</a> , <a href="#">BI985474</a>
UniProt ID:	<a href="#">Q9Z2I8</a>
Synonyms:	AF171077; AW556404; D6Wsu120e; G-SCS; GTPSCS; SCS-betaG
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)